

BTE HEARING SYSTEMS

P G4

Tech Level 3




BATTERY: 13

AMPLIFICATION: 75 dB

HOUSING COLORS

- Beige
- Gray
- Dark brown

FITTING RANGE



BASIC FEATURES

	Tech Level	3
Earhook		●
Battery compartment as on/off switch		●
Programmable push button		●
Programmable rocker switch		●
Level-dependent signal tones/melodies (can be activated/deactivated)		●
Can be used with either disposable or rechargeable battery		—
Telephone coil		●
Status LED, programmable		—
Personal color concept		—
Audiomatic power-on delay (can be activated/deactivated)		●
IP67-certified		●

SIGNAL PROCESSING

Anti-Feedback system		●
Noise Manager		
> Adaptive noise reduction		● (on/off)
> Wiener filter		● (on/off)
> Adaptive wind noise reduction, binaural		—
> Adaptive wind noise reduction		—
> Impulse suppressor		—
> Auto-situation adaption		—
> Selectronic		—
AudioTronic multi-microphone system		
> Panorama		—
> Directional static		●
> Automatic		—
> Adaptive		—
> Speech 360		—
> AudioFocus 360		—
> AudioSpot		—
Frequency and dynamics concept		
> TRC S		●
> Selective frequency compression		—
> HiFi functionality		—
> Sound Upgrade (in Wireless Streaming mode)		—
> Gain setting dependent on environment (only in Wireless Audio Streaming mode)		—
Programmable tinnitus function		—




AUTOMATIC FUNCTIONS

Comfort365		—
Acclimatic		—
Intelligent Acclimatic		—
Comformatic		—



TECHNICAL FEATURES

	Tech Level	3
Signal processing channels		12
Frequency channels		6
AGC channels		6
MPO channels		6
Hearing programs		4
> MusicSelect		—
> 2earPhone		—
> EchoClear/dereverberation		—
Data Logging		●
Wireless		
> AudioLink		—
> Binaural synchronization		—
> Wireless audio streaming		—
> CROS/BiCROS		—

ACCESSORIES | OPTIONS

Clip for battery compartment in red and blue for side recognition		●
Thin Tube with Open Tip		—
Smart Power		—
Audio streamer Smart Connect		—
Smart Remote remote control		—
Remote control Smart Key		○
Small earhook		—
AudioFix		—
Fitting Set - Thin Tube		—
Audio shoe set		○
Battery compartment safety		—
AutoPhone set		—

APPS

Smart Connect app		—
> with AudioDirSelect		—
Smart Remote app		○
> with AudioDirSelect		—

PROGRAMMING

Programming adapter 312		—
Programming adapter 13		●
Programming adapter 675		—
Directly with cable for programming interface		—

● = Standard equipment ○ = optional — = not available

More information on the individual features is available at www.audioservice.com

PG4

with earhook



MAXIMUM POWER OUTPUT

LE = 90 dB



MAXIMUM AMPLIFICATION

Amplification at LE = 50 dB

IEC 60118-7:2005¹⁾
ANSI S3.22-2009¹⁾



IEC 60118-0²⁾



TECHNICAL INFORMATION

MAXIMUM POWER OUTPUT

Peak value at 90 dB	134 dB	138 dB
1,600 Hz (RTF)	124 dB	132 dB
Mean value at high frequencies	127 dB	129 dB

MAXIMUM ACOUSTIC AMPLIFICATION

Peak value at 50 dB	70 dB	75 dB
1,600 Hz (RTF)	61 dB	67 dB
Mean value at high frequencies	64 dB	64 dB
Reference test gain	50 dB	58 dB

TECHNICAL FEATURES

Battery type	13	13
Battery life in hours	145	145
Frequency range	110 - 6,000 Hz	170 - 6,700 Hz
Battery consumption	1.40 mA	1.40 mA
Equivalent input sound pressure level of the inherent noise	24 dB	24 dB
Tinnitus Noiser, broadband	—	—
Listening coil sensitivity (1 mA/m)	93 dB	97 dB
Distortion		
500 Hz	3%	4%
800 Hz	2%	4%
1,600 Hz	1%	1%

¹⁾ Technical data measured in accordance with IEC 60118-7:2005 and ANSI S3.22-2009 at 2 ccm coupler

²⁾ Technical data measured in accordance with IEC 60118-0 at ear simulator



WARNING

Small parts present a choking hazard.
This device is not suitable for fitting to babies, small children or mentally disabled persons.



WARNING

The maximum output sound pressure level of the hearing systems can reach or exceed 132 dB SPL.
Risk of damage to the hearing of the wearer. Ensure that the hearing systems are fitted with care.

For control ranges and more programming features see Hearing System Simulation of Connex 8.4, AudioFit 8.4.5 or higher.

AS AUDIO-SERVICE GmbH · Alter Postweg 190 · 32584 Löhne · Germany info@audioservice.com · www.audioservice.com