



## SMART KEY

### THE REMOTE CONTROL THAT FITS IN YOUR POCKET

Small and simple – this is Smart Key, the new remote control. With **Smart Key**, you can switch your hearing systems on and off, change hearing programs and individually adjust volume. A key lock prevents the push buttons from being accidentally pressed when carrying the Smart Key in your pocket.

Smart Key works with high-frequency tones and is thus suitable for hearing systems with and without wireless functionality.



## SMART REMOTE

### THE CLASSIC REMOTE CONTROL

The state-of-the-art design emphasizes the quality of the **Smart Remote**. The large keys and clearly arranged color display make this remote control particularly user-friendly.

Smart Remote offers you the following functions: changing hearing programs, setting the volume, switching the hearing systems on and off, and an alarm function. A key lock prevents accidental changes to the settings of the Smart Remote.




## SMART REMOTE APP

### USE YOUR SMARTPHONE AS A REMOTE CONTROL

With our free apps, it's even simpler to operate your hearing systems.

The **Smart Remote app** uses high-frequency sounds and works with no additional accessories. Use your smartphone to control your hearing systems. This enables you to individually adjust hearing programs or the volume without anyone noticing; the latest technology even allows you to also adjust the direction and breadth of the focused area you wish to listen to.



Android is a registered trademark of Google Inc.

## EXPERIENCE INTELLIGENT TECHNOLOGY WHICH CONNECTS


Everyone's hearing is different. For this reason, **Audio Service** offers a comprehensive range of **accessory items**. Together with you, your hearing system specialist can select the perfect complement to your Audio Service hearing systems from this comprehensive product range and explain to you how to use the systems optimally.

**Just ask your hearing system specialist. He or she will be happy to advise and support you.**

The Bluetooth® word mark and the logo are registered trademarks of Bluetooth SIG Inc. Any use of this brand name by AS AUDIO-SERVICE GmbH is under license. Other trademarks and trade names remain the property of the respective holder.

# COMFORTABLE HEARING

## HEARING SYSTEM ACCESSORIES FROM AUDIO SERVICE



Art No. 105 34 864 / 10/2016 / subject to technical modifications and errors.  
Colors may vary depending on the printing process.

AS AUDIO-SERVICE GmbH · Alter Postweg 190 · 32584 Löhne · Germany  
www.audioservice.com

 **Audio Service**





## HEARING SYSTEM ACCESSORIES

### FOR EVEN GREATER COMFORT

Audio Service accessories allow you to utilize the functions of selected hearing systems to their true potential.

- **Would you like to make convenient telephone calls?**
- **Do you like to watch television?**
- **Do you enjoy listening to music?**
- **Would you like to control your hearing systems with your smartphone?**


Audio Service offers you an extensive range of accessories: remote control, charging station, audio streamer, and apps. You can even control many G4 and G3 products using your smartphone without any additional accessories.

Just ask your hearing system specialist which accessories are available for your hearing system.





## SMART CONNECT

### THE MULTI-FUNCTIONAL REMOTE CONTROL



With **Smart Connect** from Audio Service, you can use all modern infotainment options such as wireless telephone calling, listening to music, and watching television. In addition, Smart Connect also functions as a remote control.



## SMART CONNECT APP

### PERSONALIZED HEARING SYSTEM CONTROLS

Using the **Smart Connect app**, you have the option of controlling the directional characteristics of your hearing systems yourself, depending on the technical features, for example, in Tech Levels 12 and 16 the direction and breadth of the area in focus. This can be an enormous advantage in particularly challenging hearing situations.

The app is an add-on to Smart Connect. The combination of remote control and app allows you to operate your hearing systems intuitively and simply.

**Particularly practical:** If Smart Connect is not available, you can also use your smartphone as a remote control with the Smart Connect app.

Download the free  
**Smart Connect app** now



## SMART TRANSMITTER

### SUPPORTS WIRELESS CONNECTIONS

The **Smart Transmitter** is an accessory for the Smart Connect remote control so that its audio streaming function can be used with devices which do not have any Bluetooth® wireless technology.

The transmitter connects the television, stereo systems and other electronic devices which themselves have no Bluetooth function with Smart Connect. From there, the audio signals are transmitted wirelessly to your hearing systems.



## SMART POWER

### THE CHARGING STATION FROM AUDIO SERVICE

Environmental friendliness is taking on ever greater importance – we also want to offer you an environmentally-friendly alternative to standard hearing system batteries.

With **Smart Power**, constant changing of hearing system batteries is a thing of the past. Charge your hearing systems simply by using the attractive and practical charging station, avoiding the need to change the batteries constantly.

### ADVANTAGES OF SMART POWER

- + **Practical charging station**
- + **Environmentally friendly**
- + **Charges and dries your hearing systems overnight**
- + **Attractive and robust**