



# Overview In-the-ear hearing systems and accessories

October 2018



**Icon Precise<sup>1)</sup>**

**Icon**

**Sina**

**Ida BT**

Design: DIC

Design: Mini-CIC

Design: CIC

Design: Canal

Battery size 10

Battery size 10

Battery size 10

Battery size 312

With push button

With push button and AudioTronic multi-microphone system

**Equipment**

**Tech Level**

	Product	Amplification in dB	Product	Amplification in dB	Product	Amplification in dB	Product	Amplification in dB	Program change	Volume control	Wireless functionality	Remote control Smart Key	Smart Mic	Smart Transmitter 2,4	CROS quiX G5*	CROS RIC G5	Smart Direct App	Smart Remote App
<b>16</b>	Icon 16 G5 Precise	50							<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	-	-	-	-	-	<input type="checkbox"/>
			Icon 16 G5	50					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
			Icon 16 G5	55					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 16 G5	50			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 16 G5	55			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 16 G5	65			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
								Ida BT 16 G5	55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>
							Ida BT 16 G5	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>	-
<b>12</b>	Icon 12 G5 Precise	50							<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	-	-	-	-	-	<input type="checkbox"/>
			Icon 12 G5	50					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
			Icon 12 G5	55					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 12 G5	50			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 12 G5	55			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 12 G5	65			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
								Ida BT 12 G5	55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>
							Ida BT 12 G5	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>	-
<b>8</b>	Icon 8 G5 Precise	50							<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	-	-	-	-	-	<input type="checkbox"/>
			Icon 8 G5	50					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
			Icon 8 G5	55					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 8 G5	50			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 8 G5	55			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
					Sina 8 G5	65			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	-	<input type="checkbox"/>	-	-	<input type="checkbox"/>
								Ida BT 8 G5	55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>
							Ida BT 8 G5	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>	-

**Explanations**

● As standard ○ As option □ Available with optional accessories ■ Available with optional accessories, standard with Ida \* Only for in-the-ear hearing systems with wireless functionality <sup>1)</sup> Kits not available



Icon DIC	Icon Mini-CIC	Sina	Vega	Ida
Design: DIC	Design: Mini-CIC	Design: CIC	Design: Canal	Design: Canal
Battery size 10	Battery size 10	Battery size 10	Battery size 312	Battery size 312
	With push button	With push button	With push button	With rocker switch and AudioTronic multi-microphone system

**Equipment**

Tech Level	Icon DIC		Icon Mini-CIC		Sina		Vega		Ida		Program change	Volume control	Wireless functionality	Remote control Smart Key	Remote control Smart Remote*	Smart Transmitter	Smart Remote App	Smart Connect*/ Smart Connect App
	Product	Amplification in dB	Product	Amplification in dB	Product	Amplification in dB	Product	Amplification in dB	Product	Amplification in dB								
16	Icon 16 G4	40									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
	Icon 16 G4	50									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
			Icon 16 G4	40	Sina 16 G4	40	Vega 16 G4	50/55	Ida 16 G4	50/55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 16 G4	50	Sina 16 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 16 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 16 G4	65	Vega 16 G4	65	Ida 16 G4	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
12	Icon 12 G4	40									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
	Icon 12 G4	50									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
			Icon 12 G4	40	Sina 12 G4	40	Vega 12 G4	50/55	Ida 12 G4	50/55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 12 G4	50	Sina 12 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 12 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 12 G4	65	Vega 12 G4	65	Ida 12 G4	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
8	Icon 8 G4	40									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
	Icon 8 G4	50									<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	-	-	<input type="radio"/>	-
			Icon 8 G4	40	Sina 8 G4	40	Vega 8 G4	50/55	Ida 8 G4	50/55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 8 G4	50	Sina 8 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 8 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 8 G4	65	Vega 8 G4	65	Ida 8 G4	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6			Icon 6 G4	40	Sina 6 G4	40	Vega 6 G4	50/55	Ida 6 G4	50/55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 6 G4	50	Sina 6 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 6 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 6 G4	65	Vega 6 G4	65	Ida 6 G4	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 6 G4	71	Vega 6 G4	71	Ida 6 G4	71	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 6 G4	71	Vega 6 G4	71	Ida 6 G4	71	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4			Icon 4 G4	40	Sina 4 G4	40	Vega 4 G4	50/55	Ida 4 G4	50/55	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 4 G4	50	Sina 4 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 4 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 4 G4	65	Vega 4 G4	65	Ida 4 G4	65	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 4 G4	71	Vega 4 G4	71	Ida 4 G4	71	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 4 G4	71	Vega 4 G4	71	Ida 4 G4	71	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3			Icon 3 G4	40	Sina 3 G4	40	Vega 3 G4	50/55			<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
			Icon 3 G4	50	Sina 3 G4	50					<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 3 G4	55					<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 3 G4	65	Vega 3 G4	65			<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 3 G4	71	Vega 3 G4	71			<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
					Sina 3 G4	71	Vega 3 G4	71			<input checked="" type="checkbox"/>	<input type="checkbox"/>	-	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

**AS AUDIO-SERVICE GmbH**  
Alter Postweg 190  
32584 Löhne – Germany

Fon +49 5732 6878-4514  
Fax +49 5732 6878-4544  
export@audioservice.com

[www.audioservice.com](http://www.audioservice.com)  
 [facebook.com/AudioServiceGermany](https://facebook.com/AudioServiceGermany)